




Last updated on Tuesday, 09 February 2010

Vision Overview of Data Collection Screens

EScro - Enhanced Services - Lothian
 Escro - Enhanced Service - Alcohol BI - Lothian V2.00
 ESCro - Enhanced Service - Drug Misuse - Lothian V2.00

With a [consultation open](#) and [guideline tab](#) selected in Vision, click on the [EScro heading](#) to display all services from EScrio. This screen may vary depending on your Vision setup.

Double click on the guideline you want to open, and click on the button at the top of the guideline to open the EScrio data collection screens. EScrio Alcohol BI Screen




The advantage of using the data collection screens is that they provide a [bullet point summary](#) of the [enhanced service requirements](#) and coloured indicators to guide you through the consultation.

Please note: Escro data collection does not support multiple patient consultations, only one consultation should be open at any one time.

Clicking the button at the top of the guideline will launch the data collection screens. You will be prompted for your [Vision username and password](#) the [first time](#) you use the screens.

This creates a **separate consultation** and you should select the consult type and date from the dropdowns. The date can be changed if entering data retrospectively.

Consult type Surgery consultation Consult Date 08 February 2010



The summary grid provides a list of data required for this enhanced service. The grid only displays the [most recent historical entry](#) for each item.

You can navigate to the appropriate screen for data entry by clicking on the item in the grid or using the menu on the left. Any [relevant documentation or links](#) can also be found on this menu.

-  Indicates the data has never been recorded
-  Indicates the data has been recorded but is about to be or is already out of date
-  Indicates the data has been recorded and is up to date

The [coloured indicators](#) show the status of the data and the indicators on the grid represent the status on entering the consultation.



On the [data entry screens](#) the indicators will change according to the data entered during the consultation. This may be to confirm data [has been recorded](#) or to remind you that [further data is now required](#).

When the consultation is complete, you should [save the data](#) using the 'save & exit' button on the toolbar. Clicking cancel will discard any data entered.



If you need to change your vision username or password within EScro, click on settings on the toolbar. You must ensure that both the username and password are exactly the same as when logging in to Vision.